


■ LCD DISPLAY



■ WARNING INDICATOR



Battery power is low.
Please install new battery.

■ AUTO STANDBY / SHUT OFF

LCD will shut off to standby if idling for 6 seconds.
Press button to activate it.

4

This product measures body height by ultrasonic ranging technology.

■ TECHNICAL INDEX

- 1.Measuring range: 50cm - 200cm;
Measuring error: ± 1 cm; Division: 0.1cm
- 2.Power: 3x AAA batteries (excluded); Standby current: $< 20\mu A$;
Maximum operating current: 30mA.
- 3.LCD display; Buzzer reminder; One button control.

■ KEY ILLUSTRATION




All-in-one Button: One button can realize multi-functions.

- 1.Measuring: Slightly press to measure
- 2.Wake up: Slightly press to wake up from sleeping mode;
Automatic enter into sleeping mode after idling for several seconds.
- 3.Unit Conversion: Long press the button until LCD displays 9999 and release. Then the product will enter into unit conversion mode. Select the unit by slightly pressing the button; The product will automatically back out after idling for 3 seconds.

1

■ TO REPLACE BATTERY

This product use 3x1.5V AAA batteries (excluded). LCD will show Lo when battery power is low. Please install battery as described.



1. Open the battery case cover at the bottom of product
2. Take out the used batteries (with the aid of an object with pointed end if needed).
3. Install new batteries by putting one side of the battery down first and then the other side. (Please observe the polarity)

■ CARE AND MAINTENANCE

1. Do not use this product under wet, overheat or over-cool environment. (Temperature range $+5 \sim +35^{\circ}C$).
2. Clean the product with a slightly damp cloth. DO NOT immerse the product in water or use chemical or abrasive cleaning agents.
3. When in carrying, keep the product away from objects with sharp ends (such as toy, stationary etc.), in case the product is switched on frequently and consume battery power.
4. To make sure accuracy, DO NOT use this product in strong-magnetic field.
5. Keep the product away from babies and children to avoid danger of swallow the batteries and small components.
6. Take out the batteries if not used for a long time.


5

15HT-001E R0

■ HEIGHT MEASURING

Method 1: Measure height for child

1. Parents tell child to stand upright against the wall, putting the height meter over the child's head as showed in picture and holding still. Then child can walk away.
2. Keep holding the height meter and press the "●" button till hearing a beep. LCD will display the height result.



2

* The actual item may vary from the photographs/ illustrations.

CAMRY

BODY HEIGHT METER

INSTRUCTIONS FOR USE


HT-1 SERIES



ISO 9001 CERTIFIED BY SGS

Method 2: Measure height on your own

1. Keep the bottom of monitor stick closely to wall with one hand.
2. Please stand straight and put the monitor on the head flatly
3. Do not move when step back and press the button. You can see the result after hearing buzzing.



PIC 1

PIC 2

PIC 3

Remark: To make sure accuracy, there cannot be prominent barrier including skirting within the diameter range of 80cm. Please check if there is any skirting and move to the broad area when the result is not accurate. If no broad area, the real height is the result plusing height of skirting.

3