

Comparison of relevance of specific IgG₄ and specific IgG measurement in human plasma samples

The following paper summarizes the arguments of specific IgG₄ versus specific IgG testing. It should be noted that research relating to IgG₄ and its association in food-mediated diseases is still ongoing.

I Biological/Physiological Aspects

- **Regulation of IgG₄ antibodies:** As the immune response progresses there is a class switch from low-affinity antibodies to higher-affinity antibodies. The requirements for IgG₄ and IgE antibodies are similar as both require IL-4/IL-13 stimulus and both are considered part of the Th2 response. Due to the arrangement of antibody light-chains on the chromosome a B-cell can switch sequentially from a IgG₄ producing B-cell to a IgE producing B-cell (references 13-15 from publication 9).
- In contrast, IgG₁, IgG₂ and IgG₃ are less likely to play a prominent role in food-induced diseases as they are part of the Th1 immune response that is more prominent in the body's attack to antimicrobial challenges.
- Many studies report that in fact IgG₁ and IgG₃ antibodies are the only subclasses detected in response to different stimuli, which not only may prevent immediate appearance of IgG₂ and IgG₄ antibodies but also delay the class switch due to prominent IgG₁ response. This also leads to other phenomena in which IgG₄ antibodies outcompete other IgG antibodies. (1,2,3). Also, during the chronic exposure to allergens the absolute level of sIgG₄ increases and level of relative contribution of IgG₄ to total IgG is changing from 5% to up to 75% (9). This suggests that IgG₄ could be a better marker for antigens that have a long term effect on the immune system, such as foods presented in a leaky gut scenario.

- **IgG₄ and Immunotherapy:** Recent evidence has challenged the notion that successful immunotherapy is associated with elevated IgG₄ levels. For example, the Durham group showed that increased IL-10 production preceded both the clinical improvement and production of IgG₄ antibody, suggesting that IgG₄ levels are not causative for symptoms relief. The same group showed as well that the clinical benefit of immunotherapy upon treatment withdrawal is maintained, even though the allergen-specific IgG₄ levels are dropping. (4). We could conclude that these findings are not in conflict with a role of IgG₄ in mediating food intolerances.

II Pathological Aspects

There are numerous papers in various diseases suggesting a strong association or causative relation of specific IgG₄ and these pathologies.

Food elimination diets in patients with IBS/Crohns Disease: There are studies showing a positive effect for IgG₄-based food elimination diets not only for Irritative Bowel Syndrome (e.g. 8, 10, 11, 12), but also related to Crohn's disease (6,7). These papers are not summarized here as they are well known within the relevant community.

III Experimental accessibility / dynamic range of values

In an in-house study looking at the relative values of IgG and IgG₄, we tried to take literature data and our own laboratory data into account. Our study results are in line with that shown by a large study conducted in China, which suggested that food specific IgG concentrations follow more or less a normal distribution in both healthy and symptomatic adults and therefore may generate many false positive signals (5).

Identical situation can be seen in attached personal communication: in a study looking at 1674 data points measuring patient's IgG₄ and IgG specific to identical food antigens, it could be demonstrated that the dynamic range of IgG₄ is clearly superior as compared to IgG. Similar to (5), a normal distribution of IgG values can be seen, centred on Reaction Class 2 and 3,

whereas in IgG₄ values, a large body of class 0 results and a more than fourfold amount of high positives occurred. In IgG₄, distribution of classes of results do not seem to follow a normal distribution shifted to the left, but clearly distinguishable populations of negative and (highly) positive reactions are formed.

Patient specific symptoms were not correlated in this study, but taking together literature based findings and the clear indications of diagnostic opportunities in the current stories, we conclude that IgG₄ is a better marker than IgG for food intolerance testing.

References

1. Murphy SL, Li H, Mingozi F, Sabatino DE, Hui DJ, Edmonson SA, et al.; Diverse IgG subclass responses to adeno-associated virus infection and vector administration. *J Med Virol* (2009) 81(1):65–74.
2. Spinsanti LI, Farias AA, Aguilar JJ, delPilar Diaz M, Contigiani MS.; Immunoglobulin G subclasses in antibody responses to St.Louis encephalitis virus infections. *Arch Virol* (2011) 156(10):1861–4.
3. Andrew M. Collins and Katherine J. L. Jackson; A temporal model of human IgE and IgG antibody function, *Frontiers in Immunology*, 2013, Vol. 4, art. 235
4. Nouri-Aria, KT et al. Grass pollen immunotherapy induces mucosal and peripheral IL-10 responses and blocking IgG4 activity. *J Immunol* 2004; 172: 3252-9
5. Zeng Q1, Dong SY et al.; Variable food-specific IgG antibody levels in healthy and symptomatic Chinese adults; *PLoS One*, 2013, Vol. 8, Issue 1, e53612
6. N. Rajendran, D. Kumar; Food-specific IgG4-guided exclusion diets improve symptoms in Crohn's disease: a pilot study; *Colorectal Disease*, 2011, Vol. 13, Issue 9, 1009-1013
7. 3. M. Mizuno et al.; Increased incidence of allergic disorders and elevated food - specific serum IgG4 levels in Japanese patients with Crohn' s disease; *Allergology International*, 1999, Vol. 48 Issue 4, 247-251
8. Ringel-Kulka T, van Tilburg M AL, et al.; Food-specific IgG4 Antibody Titers in subjects with Food Hypersensitivity, *Gastroenterology*, 2013, Vol. 144, Issue 5, Suppl. 1, 928
9. R.C. Aalberse, S.O. Stapel et al.; Immunoglobulin G4: an odd antibody, *Clinical & Experimental Allergy*, 2009, Vol. 39, 469-477
10. Atkinson W, et al; Food elimination based on IgG antibodies in irritable bowel syndrome: a randomised controlled trial. *Gut*. 2004 Oct;53(10):1459-64.
11. Bernardi D, et al.;Time to reconsider the clinical value of immunoglobulin G4 to foods? *Clin Chem Lab Med*. 2008;46(5):687-90.
12. Zar S; Food-specific IgG4 antibody-guided exclusion diet improves symptoms and rectal compliance in irritable bowel syndrome. *Scand J Gastroenterol*. 2005 Jul;40(7):800-7.