





		Distribution of Age							
		3~10	11~18	19~29	30~39	40~49	50~59	60~	Total
MALE	N	4	11	33	24	19	15	9	115
FEMALE	N	3	18	30	14	29	21	9	124
	Total	7	29	63	38	48	36	18	239

		Distribution of Age							
		3~10	11~18	19~29	30~39	40~49	50~59	60~	
	N	4	11	33	24	19	15	9	
MALE	Height(cm)	122.8±14.5	170.5±9.4	176.4±6.5	173.5±4.9	172.8±5.7	167.5±7.1	167.5±4.7	
	Weight(kg)	28.6±12.8	72.6±20.4	74.9±12	78.6±14.1	73.5±9.7	69±4.8	69.2±10.3	
	BMI(kg/m <sup>2</sup> )	17.8±3.6	24.8±5.8	24.1±3.7	26±4.1	24.6±3.1	24.7±2.1	24.6±2.5	
	N	3	18	30	14	29	21	9	
FEMALE	Height(cm)	128.9±12.3	161.3±4.3	159.4±4.3	159.3±5.1	159±5.5	155.7±4.7	154.7±5.1	
	Weight(kg)	30±6.4	57.4±9.9	58.1±8.3	58.4±9.1	59.6±9.5	58.7±6.9	58.5±6	
	BMI(kg/m <sup>2</sup> )	17.8±0.4	22.1±3.8	22.9±3.3	23±3.4	23.6±3.5	24.2±2.8	24.5±2.4	

Section	R <sup>2</sup> (SEE)	Correlation graph
Fat-free Mass	0.984 (1.485)	 <p>Simple Scatter with Fit Line of FFM_Total_DEXA by FFM_Total_I_Series</p> <p>R<sup>2</sup> Linear = 0.984 SEE = 1.485kg</p>
Muscle Mass	0.983 (1.391)	 <p>Simple Scatter with Fit Line of Muscle_DEXA by Muscle_I_Series</p> <p>R<sup>2</sup> Linear = 0.983 SEE = 1.391kg</p>



## 8. References

1. Thomasset MA., Bioelectric properties of tissue. Impedance measurement in clinical medicine. Significance of curves obtained. Lyon Med. 1962 Jul 15;94:107-18.
2. Earl C. Hoffer et al., Correlation of whole-body impedance with total water volume. Journal of applied physiology Vol.27, No.4, October 1969
3. CLINICAL EVALUATION: A GUIDE FOR MANUFACTURERS AND NOTIFIED BODIES UNDER DIRECTIVES 93/42/EEC and 90/385/EEC, MEDDEV 2.7/1 revision 4, June 2016
4. World Medical Association Declaration of Helsinki: Ethical principles for medical research involving human subjects
5. Despres, J. P., Moorjani, S., Lupien, P. J., Tremblay, A., Nadeau, A., & Bouchard, C. (1990). Regional distribution of body fat, plasma lipoproteins, and cardiovascular disease. *Arteriosclerosis, Thrombosis, and Vascular Biology*, 10(4), 497-511.
6. Siri WE., Body composition from fluid spaces and density\_Analysis of methods, *Nutrition*. 1993 Sep-Oct;9(5):480-91; discussion 480, 492.
7. Wang ZM, Pierson RN Jr, Heymsfield SB., The five-model a new approach to organizing body-composition research, *Am J Clin Nutr*. 1992 Jul;56(1):19-28.
8. Vivian H. Heyward, Dale R. Wagner, *Applied Body Composition Assessment - 2nd*. 2004 Human Kinetics
9. Yoshizumi T, Nakamura T, Yamane M, Islam AH, Menju M, Yamasaki K, Arai T, Kotani K, Funahashi T, Yamashita S, Matsuzawa Y., Abdominal fat: Standardized Technique for Measurement at CT. *Radiology*. 1999 Apr;211(1):283-6.
10. Lukaski HC, Johnson PE, Bolonchuk WW, Lykken GI., Assessment of fat-free mass using bioelectrical impedance measurements of the human body. *Am J Clin Nutr*. 1985 Apr;41(4):810-7
11. Hansen RD, Raja C, Aslani A, Smith RC, Allen BJ., Determination of skeletal muscle and fat-free mass by nuclear and dual-energy X-ray absorptiometry methods in men and women aged 51-84 y. *Am J Clin Nutr*. 1999 Aug;70(2):228-33.
12. Atkin LM, Davies PS., Diet composition and body composition in preschool Children. *Am J Clin Nutr*. 2000 Jul;72(1):15-21.
13. J. C. MorkebergH.-P. ShengW. W. Wong, Extracellular Water Estimated by the bromide dilution method from samples of urine, saliva, and plasma. In *Vivo Body Composition Studies* pp 55-59